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Department:		Room:	
Material-No.:		Serial-No.:	
Contract-No.:		Expire date:	
Order-No.:		System - ID:	

The instructions C2-015.831.01.10.XX are required for this protocol

Evaluating the Condition of the System

The system has no deficiencies. *	☐
The system has slight deficiencies that have no effect on continued operation of the system. The deficiencies should be corrected preventively. *	☐
The system has serious deficiencies. For safety reasons, continued operation of the system is permitted only after successfully correcting the deficiencies.	☐

*) For imaging systems: The result of the image quality check shows no deviation from the reference values.

The evaluation was performed after completing all work steps.

Date: _____ Name: _____

If required by country-specific regulations:

The customer or a representative has taken note of the result of the evaluation of the system condition.

The customer or a representative has taken note of the result of the evaluation of the system condition.

Remarks:

Remarks Regarding the Protocol

The chapter numbers refer to the chapters in the instruction, which is referenced on the cover page.

All pages have to have the serial number of the system and the date of maintenance in the page header.

The assignment n.a. (not applicable) indicates that the checkpoint or measured value is not used for this system.

On page 2 the completeness and the results of the maintenance work is confirmed.

Explanation of Abbreviations in the Protocol

Abbrev.	Explanation	Abbrev.	Explanation
SI	Safety Inspection	PMF	Preventive Maintenance, Operating Value Check, Function Check
SIE	Electrical Safety Inspection	Q	System Quality, Image Quality
SIM	Mechanical Safety Inspection	QIQ	Image Quality
PM	Preventive Maintenance	QSQ	System Quality Check
PMP	Periodic Preventive Maintenance	SW	Software Maintenance
PMA	Preventive Maintenance Adjustments	CSE	Customer Service Engineer

Measuring Equipment and Measuring Instruments Used

Measuring instruments and measuring devices (phantoms, MR coils, etc.) may not be entered in the table if they have already been entered in the mobile device.

Measuring equipment / -instruments	Serial No.	Date used

1 General

1.9 Special information for Mobile CT & Mobile Stroke Unit

2 Partial Maintenance (every 6 months)

2.1 Switching on the System

2.2 Switching off the System

2.3 Gantry: Removing the covers

2.4 PHS1A: Cover Removal & Checking the Horizontal Drive Toothed Belt Tension

PM Horizontal drive toothed belt tension checked (6M) ==>

☐☐☐

NOTE: This check should be performed only once during the first maintenance after system/PHS installation or after horizontal drive or toothed belt replacement.

Value:

2.5 Mobile CT/Mobile Stroke Unit (optional): Inspections for Mobile CT/Mobile Stroke Unit

SI Mobile kit inspected for damage (6M)

☐☐☐

SI Patient table inspected for damage and loose parts (6M)

☐☐☐

SI Gantry inspected for damage, loose parts and screw connections (6M)

☐☐☐

2.6 Sliding gantry (option): Checking the switching elements

PM Switching elements at front and rear checked (6M)

☐☐☐

2.7 Gantry: Cleaning the slip ring assembly

PM Slip ring assembly cleaned (6M)

☐☐☐

2.8 Gantry: Checking/cleaning/replacing the Power Carbon Brushes (565VDC)

PM Power carbon brushes (565V DC) checked (6M)

☐☐☐

PM Power carbon brushes (565V DC) replaced (6M)

☐☐☐

2.9 Gantry: Checking/cleaning/replacing the Power Carbon Brushes (230VAC)

PM Power carbon brushes (230V AC) checked (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

PM Power carbon brushes (230V AC) replaced (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

OK not n.a.
 O
 K

2.10 Gantry: Checking/cleaning/replacing the Signal Brush Block

PM Signal brush block checked & cleaned (6M)

☐	☐	☐
☐	☐	☐

PM Signal brush block replaced (6M)

2.11 Gantry: Checking/replacing the MAS Air Filter

PM MAS air filter checked (6M)

☐	☐	☐
☐	☐	☐

PM MAS air filter replaced or washed (6M)

2.12 IMS: Cleaning the air inlets of the ICS and IRS

PM Air inlets of ICS and IRS cleaned (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

2.13 System: Reinstalling the covers

SIE Covers and protective ground wires installed (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

2.14 Switching on the System

2.15 System: Cleaning the System

PM System cleaned (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

2.16 System: Quality assurance test

Q Quality assurance test completed (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

2.17 System: Checking the system functionality

Q Functionality of the CT system checked (6M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

2.18 Maintenance Documentation

Work steps of this section were performed by:

Signature: _____

Date: _____

Name: _____

3 Full Maintenance (every 12 months)

3.1 Switching on the System

3.2 Gantry: Tilting the Gantry

PM Gantry tilt over the full $\pm$ tilt range executed (12M)

☐ ☐ ☐

OK not n.a.
 O
 K

3.3 Switching off the System

3.4 Gantry: Removing the covers

3.5 PHS1A: Cover Removal

3.6 System: Visual Inspection according to DIN EN 62353/IEC 62353

SI	Gantry: Visually inspected (12M)	☐	☐	☐
SI	Patient table: Visually inspected (12M)	☐	☐	☐
SI	System cabinets: Visually inspected (12M)	☐	☐	☐
SI	Balancing weights: Visually inspected (12M)	☐	☐	☐
SI	System cables & cable routing: Visually inspected (12M)	☐	☐	☐
SI	Accessories: Visually inspected (12M)	☐	☐	☐
SI	Phantoms & phantom holder: Visually inspected (12M)	☐	☐	☐
SI	Care Vision CT (option): Visually inspected (12M)	☐	☐	☐
SI	Care Vision CT cables & cable routing (option): Visually in-spected (12M)	☐	☐	☐
SI	Sliding gantry (option): Visually inspected	☐	☐	
	☐			

3.7 PHS: Measuring the leakage current

SI	Leakage current (tabletop plate to GND terminal) < 0.1 mA [AC] (12M)	☐
	Value:	
SI	Leakage current (tabletop plate to GND terminal) < 0.01 mA[DC] (12M)	☐
	Value:	
SI	Leakage current (ECG leads (option) to GND terminal) < 0.1mA [AC] (12M)	☐
	Value:	

0.01 mA [DC] (12M)

☐☐☐

Value:

SI Leakage current (ECG leads (option) to tabletop plate) <
0.1mA [AC] (12M)

☐☐☐

Value:

SI Leakage current (ECG leads (option) to tabletop plate) <

☐☐☐

0.01 mA [DC] (12M)

Value:

OK not n.a.
O
K

3.8 Mobile CT/Mobile Stroke Unit (option): Inspections for Mobile CT/ Mobile Stroke Unit

- | | | | | |
|----|---|--------------------------|--------------------------|--------------------------|
| SI | Mobile kit inspected for damage (6M) | ☐ | ☐ | ☐ |
| SI | Patient table inspected for damage and loose parts (6M) | ☐ | ☐ | ☐ |
| SI | Gantry inspected for damage, loose parts and screw connections (6M) | ☐ | ☐ | ☐ |

3.9 System: Inspecting the labels

- | | | | |
|----|-------------------------------|--------------------------|--------------------------|
| SI | System labels inspected (12M) | ☐ | ☐ |
| | ☐ | | |

3.10 PHS1A: Lubricating the Vertical Drive Spindle

- | | | | |
|----|---|--------------------------|--------------------------|
| PM | Vertical drive spindle lubricated (12M) | ☐ | ☐ |
| | ☐ | | |

3.11 PHS1A: Lubricating the Guide Rails of the Table

- | | | | |
|----|------------------------------|--------------------------|--------------------------|
| PM | Guide rails lubricated (12M) | ☐ | ☐ |
| | ☐ | | |

3.12 PHS1A: Checking the Horizontal Drive Toothed Belt Tension

- | | | | |
|----|--|--------------------------|--------------------------|
| PM | Horizontal drive toothed belt tension check (6M) ==> NOTE: | ☐ | ☐ |
| | ☐ | | |

This check should be performed only once during the first
maintenance after system/PHS installation or after horizontal drive
or toothed belt replacement.

Value:

3.13 PHS1A: Checking the Anchoring of the PHS

- | | | | | |
|----|-------------------------|--------------------------|--------------------------|--------------------------|
| PM | Anchoring checked (12M) | ☐ | ☐ | ☐ |
|----|-------------------------|--------------------------|--------------------------|--------------------------|

3.14 MCS (option): Checking the monitor ceiling system

- SI MCS visually checked for cracks (12M)
- PM MCS movement functions checked (12M)
- PM MCS mechanical end stops checked (12M)
- PM MCS load balance checked (12M)

3.15 Sliding gantry (option): Function Check

- PM Toothed belt checked (12M)
- PM Switching elements, front and rear checked (6M)
- PM Protective conductor connection checked (12M)

3.16 Gantry: Checking the rotation drive belt and belt tension

- PM Belt checked (12M)
- PM Belt tension checked (12M)

OK	not	n.a.
	O	
	K	

3.17 Gantry: Checking the Rotation Brake**PM** Rotation brake checked (12M)**3.18 Gantry: Checking the Cover Hooks****PM** Elasticity of the hook holders on both the front and rear of the gantry checked (12M)**3.19 Gantry: Cleaning the Air Outlets on the Gantry Top****PM** Gantry air outlets cleaned (12M)**3.20 Gantry: Checking the Gantry Cooling Fans****PM** Gantry cooling fans checked (12M)**3.21 Gantry: Cleaning the Slip Rings****PM** Slip ring assembly cleaned (6M)**3.22 Gantry: Checking/cleaning/replacing the Power Carbon Brushes (565VDC)****PM** Power carbon brushes (565V DC) checked (6M)**PM** Power carbon brushes (565V DC) replaced (6M)**3.23 Gantry: Checking/cleaning/replacing the Power Carbon Brushes (230VAC)****PM** Power carbon brushes (230V AC) checked (6M)**PM** Power carbon brushes (230V AC) replaced (6M)**3.24 Gantry: Checking/cleaning/replacing the Signal Brush Block****PM** Signal brush block checked & cleaned (6M)**PM** Signal brush block replaced (6M)

3.25 Gantry: Checking/replacing the Air Filters for the Gantry

PM Gantry air filters checked (12M)

☐	☐	☐
☐	☐	☐

PM Gantry air filters replaced (12M)

3.26 Gantry: Lubricating the Main Bearing

PM Main bearing lubricated (12M)

☐	☐	☐
--------------------------	--------------------------	--------------------------

3.27 Gantry: Checking/replacing the MAS Air Filter

PM MAS air filter checked (6M)

☐	☐	☐
☐	☐	☐

PM MAS air filter replaced or washed (6M)